

Instrument: BD FACSCalibur

Laser	Detector	Filter	Mirror	Some common fluorochromes
Blue - 488 nm	FL1	530/30	560SP	AF488, FITC, GFP, Zombie Green, BB515
	FL2	585/42	640LP	PE
	FL3	—	670LP	PE-Texas Red, PE-Cy5, PE-Cy5.5, PE-Cy7, PerCP, PerCP-Cy5.5, PI, 7-AAD, BB700
		488/10	—	SSC, FSC
Red - 635 nm	FL4	661/16	640LP	APC, APC-Cy7, AF700